

## Bacteria-free check media for marine algae

### MM 23 (M. Tatewaki, pers. comm.)

NaCl	1.8 g
MgSO <sub>4</sub> · 7H <sub>2</sub> O	500 mg
KCl	60 mg
NaNO <sub>3</sub>	100 mg
CaCl <sub>2</sub> · 2H <sub>2</sub> O	36.7 mg
K <sub>2</sub> HPO <sub>4</sub>	6 mg
Sucrose	400 mg
PII metals	2 mL
FeCl <sub>3</sub> · 6H <sub>2</sub> O	48 µg
Thiamine HCl	10 µg
Biotin	0.1 µg
Vitamin B <sub>12</sub>	0.2 µg
C-Source Mix II	1 mL
Tris (hydroxymethyl) aminomethane	100 mg
Distilled water	97 mL
pH 8.0	

### P II metals

Na <sub>2</sub> EDTA · 2H <sub>2</sub> O	100 mg
H <sub>3</sub> BO <sub>3</sub>	114 mg
FeCl <sub>3</sub> · 6H <sub>2</sub> O	4.9 mg
MnSO <sub>4</sub> · 4H <sub>2</sub> O	16.4 mg
ZnSO <sub>4</sub> · 7H <sub>2</sub> O	2.2 mg
CoSO <sub>4</sub> · 7H <sub>2</sub> O	480 µg
Distilled water	100 mL

### Reference

Provasoli, L. 1963 Growing marine seaweeds. In *Proceedings of the Fourth International Seaweed Symposium*, University of Tokyo Press, Tokyo, p. 9-17.

## Bacteria-free check media for marine algae

### **C-Source Mix II** (M. Tatewaki, pers. comm.)

Glycine	100	mg
D,L-Alanine	100	mg
L-Asparagine	100	mg
Sodium acetate · 3H <sub>2</sub> O <sup>1)</sup>	200	mg
Glucose	200	mg
L-Glutamic acid	200	mg
Distilled water	100	mL

1) In the NIES-Collection, 200 mg sodium acetate · 3H<sub>2</sub>O is replaced by 120 mg sodium acetate, anhydrous.