

Media for protozoa

Waris EDTA

EDTA	0.522	g
FeSO ₄ · 7H ₂ O	0.498	g
KOH (1 mol/L solution)	5.4	mL
Distilled water	94.6	mL

EDTA and FeSO₄ · 7H₂O is heated for 30 min (100 °C); KOH is added to the cooled mixture.

Reference

McFadden, G. I., Melkonian, M. 1986 Use of Hepes buffer for microalgal culture media and fixation for electron microscopy. *Phycologia*, **25**, 551-557.