mTYGM-9

Pre-solution

K ₂ HPO ₄	0.28 g
KH ₂ PO ₄	40 mg
Casein Digest	0.2 g
Yeast extract	0.1 g
NaCl	0.75 g
Mucin, gastric	0.2 g
Distilled water	97 mL

Sterilize Pre-solution by autoclaving (121°C, 15 min) , add aseptically 3 mL horse serum and 50 μ L 10% Tween 80 dissolved in absolute ethanol (both filter-sterilized). Keep in a cool place.