

Waris P-2

Na ₂ EDTA · 2H ₂ O	0.3	g
H ₃ BO ₃	0.114	g
MnCl ₂ · 4H ₂ O	14.4	mg
ZnSO ₄ · 7H ₂ O	2.1	mg
CoCl ₂ · 6H ₂ O	0.4	mg
Distilled water	100	mL

Reference

McFadden, G. I., Melkonian, M. 1986 Use of Hepes buffer for microalgal culture media and fixation for electron microscopy. *Phycologia*, **25**, 551-557.