Media for protozoa

BLEY (Barley leaf & Egg-yolk extracted liquid)

BL solution: 1 pack (3 g) of Barley leaf (Trade name: Oomugi Wakaba Funmatu 100%, Yamamoto Kanpou Seiyaku Co., Ltd.) is added to 1,000 mL distilled water and boiled for 10 minutes, while stirring. The solution is passed through the qualitative filter paper.

EY solution: 9 g of crushed yolk of hard-boiled egg is added to 1,000 mL distilled water and boiled for 10 minutes, while stirring. The solution is passed through the qualitative filter paper.

BL and EY Solution are mixed. The mixed solution is dilution with distilled water up to 6,000 mL, and then sterilized by autoclaving.