

IMK/2

Into 100 mL seawater dissolve 12.6 mg powder medium of Daigo IMK (Nihon Pharmaceutical Co., Ltd.).

In the NIES-Collection, IMK/2 medium is used after autoclaving (121°C, 20 min).

Reference

Takahashi, K., Lum, W. M., Benico, G., Uchida, H., Ozawa, M., Oikawa, H., Suzuki, T., Nguyen, N. V., Ha, D. V., Iwataki, M. 2021 Harmful algal bloomsToxigenic strains of *Azadinium poporum* (Amphidomataceae, Dinophyceae) from Japan and Vietnam, with first reports of *A. poporum* (ribotype A) and *A. trinitatum* in Asian Pacific. J. Phycol., 69, 175-187.