

## Media for protozoa

### PYG

PPYG medium is diluted half strength.

Proteose Peptone	0.15 g
Yeast extract	0.1 g
Glucose	0.5 g
Distilled water	100 ml
pH	7.0

In the NIES-Collection, the amount of Proteose peptone is reduced from 0.2g to 0.15g.

#### Reference

Mori, K., Yomo, T., Kashiwagi, A. 2011 Single-cell isolation and cloning of *Tetrahymena thermophila* cells with a fluorescence-activated cell sorter. *J. Eukaryotic Microbiol.*, **58**(1), 37-42.