

## Media for freshwater, terrestrial, hot spring and salt water algae

### O

Glucose	100 mg
Tryptone	100 mg
Yeast extract	100 mg
Beef extract <sup>1)</sup>	50 mg
Agar	150 mg
Distilled water	100 ml

1) In the NIES-Collection, beef extract is removed.  
Indicated as "Ochomonas medium" in reference.

#### References

Ichimura, T., Itoh, T. 1977 17. Preservation methods of microalgae (I) [17. Bisaisôruï no hozonhô (I)]. In *Preservation methods of microorganisms [Biseibutsu Hozonhô]*, Ed. by Nei, T., University of Tokyo Press, Tokyo, p. 355-373 (in Japanese without English title).  
Starr, R. C. 1964 The culture collection of algae at Indiana University. *Amer. J. Bot.*, **51**, 1013-1044.