Tre

To 100 mL MBM medium add 1 g proteose peptone and 2 g glucose.

Indicated as "Trebouxia medium" in reference.

Reference

Ichimura, T., Itoh, T. 1977 17. Preservation methods of microalgae (I) [17. Bisaisôrui no hozonhô (I)]. In *Preservation methods of microorganisms [Biseibutsu Hozonhô]*, Ed. by Nei, T., University of Tokyo Press, Tokyo, p. 355-373 (in Japanese without English title).

Starr, R. C. 1964 The culture collection of algae at Indiana University. Amer. J. Bot., 51, 1013-1044.

MBM

KNO ₃	25	mg
MgSO ₄ · 7H ₂ O	7.5	mg
K ₂ HPO ₄	7.5	mg
KH ₂ PO ₄	17.5	mg
NaCl	2.5	mg
CaCl ₂ · 2H ₂ O	1	mg
Fe solution	0.1	mL
A5 solution	0.1	mL
Agar	1.5	g
Distilled water	99.8	mL
pH 6.0		

Indicated as "Modified Bristol medium" in reference.

Reference

Ichimura, T., Itoh, T. 1977 17. Preservation methods of microalgae (I) [17. Bisaisôrui no hozonhô (I)]. In *Preservation methods of microorganisms [Biseibutsu Hozonhô]*, Ed. by Nei, T., University of Tokyo Press, Tokyo, p. 355-373 (in Japanese without English title).

Fe solution

FeSO ₄ · 7H ₂ O	200	mg
Distilled water	100	mL
Conc · H ₂ SO ₄	0.026 mL ¹⁾	

1) 2 drops/500 mL

Reference

Ichimura, T., Itoh, T. 1977 17. Preservation methods of microalgae (I) [17. Bisaisôrui no hozonhô (I)]. In *Preservation methods of microorganisms [Biseibutsu Hozonhô]*, Ed. by Nei, T., University of Tokyo Press, Tokyo, p. 355-373 (in Japanese without English title).

Media for freshwater, terrestrial, hot spring and salt water algae

A₅ solution

H ₃ BO ₃	286	mg
$MnSO_4 \cdot 7H_2O^{1)}$	250	mg
ZnSO ₄ ·7H ₂ O	22.2	mg
CuSO ₄ · 5H ₂ O	7.9	mg
Na ₂ MoO ₄ · 2H ₂ O	2.1	mg
Distilled water	100	mL

¹⁾ In the NIES-Collection, 250 mg MnSO $_4 \cdot$ 7H $_2$ O is replaced by 217 mg MnSO $_4 \cdot$ 5H $_2$ O.

Reference

Holm-Hansen, O., Gerloff, G. C., Skoog, F. 1954 Cobalt as an essential element for blue-green algae. *Physiol. Planta.*, **7**, 665-675.